


Product Datasheet

AB-MECA (orb341594)

Description	AB-MECA is a high affinity A3 adenosine receptor agonist, has high affinity for recombinant A1 and A...
Note	For research use only
Formula	C ₁₈ H ₂₁ N ₇ O ₄
MW	399.403842687607
CAS Number	[152918-26-8]
SMILES	<chem>NC1C=CC(CNC2=NC=NC3N([C@@H]4O[C@H](C(NC)=O)[C@@H](O)[C@H]4O)C=NC2=3)=CC=1</chem>
Expiration Date	12 months from date of receipt.



Chemical
structure of
AB-MECA