


Product Datasheet

Sal003 (orb341656)

Description	Sal003 is a potent cell-permeable analog of the eIF2α phosphatase inhibitor Salubrinal with enhance...
Note	For research use only
Formula	C ₁₈ H ₁₅ Cl ₄ N ₃ OS
MW	463.21
CAS Number	[1164470-53-4]
SMILES	O=C(C=CC1C=CC=CC=1)NC(C(Cl)(Cl)Cl)NC(=S)NC1C=CC(Cl)=CC=1
Expiration Date	12 months from date of receipt.



Chemical
structure of
Sal003