

## Product Datasheet

### Pimavanserin(ACP-103) (orb341711)

|                        |                                                                                                                                 |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Pimavanserin(ACP-103) is a potent and selective 5-HT <sub>2A</sub> receptor inverse agonist with mean pIC <sub>50</sub> of w... |
| <b>Note</b>            | For research use only                                                                                                           |
| <b>Formula</b>         | C <sub>25</sub> H <sub>34</sub> FN <sub>3</sub> O <sub>2</sub>                                                                  |
| <b>MW</b>              | 427.56                                                                                                                          |
| <b>CAS Number</b>      | [706779-91-1]                                                                                                                   |
| <b>SMILES</b>          | <chem>N(CC1=CC=C(F)C=C1)(C1CCN(C)CC1)C(NCC1=CC=C(OCC(C)C)C=C1)=O</chem>                                                         |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                 |