

Product Datasheet

SirReal2 (orb341780)

Description	SirReal2 is a potent and selective Sirt2 inhibitor with IC50 of 140 nM.
Note	For research use only
Formula	C ₂₂ H ₂₀ N ₄ OS ₂
MW	420.55
CAS Number	[709002-46-0]
SMILES	<chem>CC1=CC(C)=NC(SCC(NC2=NC=C(CC3=CC=CC4C=CC=CC3=4)S2)=O)=N1</chem>
Expiration Date	12 months from date of receipt.