


Product Datasheet

TG003 (orb341894)

Description	TG003 is a potent and ATP-competitive Cdc2-like kinase (Clk) inhibitor with IC50 of 20 nM, 200 nM, a...
Note	For research use only
Formula	C ₁₃ H ₁₅ NO ₂ S
MW	249.3287
CAS Number	[300801-52-9]
SMILES	<chem>COC1C=C2N(CC)\C(=C\C(C)=O)SC2=CC=1</chem>
Expiration Date	12 months from date of receipt.



Chemical
structure of
TG003