

Product Datasheet

Paox antibody (FITC) (orb34206)

Description	Paox antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Peroxisomal N(1)-acetyl-spermine/spermidine oxidase protein (1-504AA)
Target	Paox
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8C0L6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.