

Product Datasheet

Epitope Tag Antibody Sampler Kit (orb342265)

Description

Epitope Tag Antibody Sampler Kit

Conjugation

Unconjugated

Tested

ChIP, ELISA, FC, IF, IHC, IP, WB

Applications
Immunogen

GFP, GST, DYKDDDDK peptide (the same epitope recognized by Sigma's Anti-FLAG antibodies), 6X His peptide

Storage

Store Epitope Tag Kit at -20° C prior to opening. Aliquot antibodies and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Epitope Tag kit contains: 6X HIS Antibody 100 µg DYKDDDDK Antibody (Binds to same epitope as Sigma's Anti-FLAG M2 Antibody) 250 µg GFP Antibody 100 µg GST Antibody 1.0mg Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated 100 µg Anti-GOAT IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated 100 µg

Purity

6X HIS antibody, GFP Antibody, DYKDDDDK Antibody (Binds to same epitope as Sigma's Anti-FLAG M2 Antibody), GST Antibody and Conjugated Secondary Antibodies are all affinity purified.

Dilution Range

ELISA: User Optimized, ChIP: User Optimized, FC: User Optimized, IHC: User Optimized, IF: User Optimized, IP: User Optimized, WB: User Optimized

Expiration Date

Please enquire.