

Product Datasheet

Western Blot Kit (orb342271)

Description	Western Blot Kit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!
Form/Appearance	n/a
Concentration	n/a
Storage	See kit insert for complete instructions.
Note	For research use only
Application notes	Biorbyt' Chemiluminescent Western Blot Kit, namely our FemtoMax Chemiluminescent Substrate Kit for Western Blotting, combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol substrate. The Chemiluminescent Western Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains sufficient substrate for up to 30 mini blots at 7.5 x 8 cm ² and is stable for at least 1 year when stored at +4°C.
Purity	This kit contains: anti-mouse peroxidase, mouse IgG control, BSA, tween-20, and Femtomax, along with the protocol and additional supplies. See kit insert for complete details and instructions.
Dilution Range	WB: 1:20,000
Expiration Date	Please enquire.