

Product Datasheet

Biotin Anti-Rabbit IgG Antibody Maxtag Histo For IHC (orb342290)

Description	Biotin Anti-Rabbit IgG Antibody Maxtag Histo For IHC
Conjugation	Biotin
Tested Applications	IHC
Immunogen	Rabbit IgG whole molecule
Preservatives	0.01% (w/v) Sodium Azide and 0.01% (w/v) Gentamicin Sulfate
Form/Appearance	Liquid (sterile filtered)
Concentration	1.0 mg/mL
Storage	See kit insert for complete instructions.
Note	For research use only
Application notes	See Kit insert for complete instructions.
Expiration Date	12 months from date of receipt.