

## Product Datasheet

# Biotin anti-Goat IgG Antibody MaxTag Histo for IHC (orb342291)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Biotin anti-Goat IgG Antibody<br>MaxTag Histo for IHC          |
| <b>Conjugation</b>         | Biotin                                                         |
| <b>Tested Applications</b> | IHC                                                            |
| <b>Immunogen</b>           | Goat IgG whole molecule                                        |
| <b>Preservatives</b>       | 0.01% (w/v) Sodium Azide and<br>0.01% (w/v) Gentamicin Sulfate |
| <b>Form/Appearance</b>     | Liquid (sterile filtered)                                      |
| <b>Storage</b>             | See kit insert for complete<br>instructions.                   |
| <b>Note</b>                | For research use only                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                |