

Product Datasheet

MaxTag Histo for IHC (orb342293)

Description	Formaldehyde Fixative 16% (w/v) Ultrapure EM Grade (5 pack - 10 ml per ampule) MaxTag Histo for IHC
Conjugation	Unconjugated
Tested Applications	IHC
Form/Appearance	Liquid
Storage	See kit insert for complete instructions.
Note	For research use only
Application notes	Formaldehyde Fixative 16% (w/v) Ultrapure EM Grade (5 pack - 10 ml per ampule) MaxTag Histo for IHC is designed for use in immunohistochemistry experiments.
Dilution Range	IHC: User Optimized
Expiration Date	12 months from date of receipt.