

Product Datasheet

Polymount Mounting Media Maxtag Histo For IHC (orb342295)

Description	Polymount Mounting Media Maxtag Histo For IHC
Conjugation	Unconjugated
Tested Applications	IF, IHC
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid
Concentration	Neat
Storage	See kit insert for complete instructions.
Note	For research use only
Application notes	This component is ready-to-use and requires no additional preparation.
Dilution Range	IHC: See kit insert., IF: See kit insert.
Expiration Date	12 months from date of receipt.