

Product Datasheet

Rat Peroxidase Anti-Peroxidase (PAP) Soluble Complex (orb342322)

Description

Rat Peroxidase Anti-Peroxidase (PAP)

Conjugation

HRP

Tested

IHC

Applications
Immunogen

Peroxidase (Horseradish)

Preservatives

Preservative: None. Stabilizer: None.
0.01 M Sodium Phosphate, 0.15 M
Sodium Chloride, pH 7.2

Form/Appearance

Lyophilized

Concentration

5.0 mg/mL

Storage

Store vial at 4° C prior to restoration.
For extended storage aliquot
contents and freeze at -20° C or
below. Avoid cycles of freezing and
thawing. Centrifuge product if not
completely clear after standing at
room temperature. This product is
stable for several weeks at 4° C as
an undiluted liquid. Dilute only prior
to immediate use.

Note

For research use only

Application notes

Suitable for use in peroxidase
enzymatic labeling experiments.
After reconstitution, dilute sample
with PBS containing 5% normal
serum from same host bridging
antibody.

Purity

This product was prepared from rat
ascites by delipidation, salt
fractionation and ion exchange
chromatography. Assay by
immunoelectrophoresis resulted in a
single precipitin arc against anti-Rat
Serum and anti-Peroxidase.

Dilution Range

IHC: User Optimized

Expiration Date

12 months from date of receipt.