

Product Datasheet

Sheep Peroxidase Anti-Peroxidase (PAP) Soluble Complex (orb342324)

Description

Sheep Peroxidase Anti-Peroxidase

Conjugation	HRP
Tested Applications	IHC
Immunogen	Peroxidase (Horseradish)
Preservatives	Preservative: None. Stabilizer: None. 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form/Appearance	Lyophilized
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Suitable for use in peroxidase enzymatic labeling experiments. After reconstitution, dilute sample with PBS containing 5% normal serum from same host bridging antibody.
Purity	This product was prepared from Sheep antiserum by delipidation, salt fractionation and ion exchange chromatography. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Sheep Serum and anti-Peroxidase.
Dilution Range	IHC: User Optimized
Expiration Date	12 months from date of receipt.