

Product Datasheet

NFkB Antibody Sampler Kit (orb342330)

Description	NFkB Antibody Sampler Kit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	Liquid (sterile filtered)
Storage	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Each kit contains one unit of the following products: Anti-NFkB p50 (NFkB1) (RABBIT) Antibody (100 µL) Anti-NFkB p65 (Rel A) (RABBIT) Antibody (100 µL) Anti-NFkB cRel (RABBIT) Antibody (100 µL) Anti-IkBα C-terminal (RABBIT) Antibody (100 µL) Anti-RABBIT IgG (H&L) (GOAT) Peroxidase Conjugated Antibody (100 µg) CONTROL PEPTIDE for Anti-NFkB p50 (NFkB1) (RABBIT) Antibody (50 µg) CONTROL PEPTIDE for Anti-NFkB p65 (Rel A) (RABBIT) Antibody (50 µg) CONTROL PEPTIDE for Anti-NFkB cRel (RABBIT) Antibody (50 µg) CONTROL PEPTIDE for Anti-IkBα C-terminal (RABBIT) Antibody (50 µg)
Dilution Range	ELISA: 1:5,000 - 1:25,000, WB: 1:500 - 1:2,000
Expiration Date	Please enquire.