

Product Datasheet

Tumor Suppressor Protein Antibody Sampler Kit (orb342333)

Description	Tumor Suppressor Protein Antibody Sampler Kit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Form/Appearance	Liquid (sterile filtered)
Storage	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Each kit contains one unit of the following products: (100 µL) p/n 100-401-173 Anti-p130 (Rb2) (RABBIT) Antibody (100 µL) p/n 100-401-170 Anti-p16 (RABBIT) Antibody (100 µL) p/n 100-401-172 Anti-p27 (RABBIT) Antibody (100 µg) p/n 611-1322-0100 Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated
Dilution Range	ELISA: 1:2,000 - 1:10,000, IP: 1:100, WB: 1:500 - 1:5,000
Expiration Date	Please enquire.