


Product Datasheet

ADGRA3 antibody (orb34237)

Description	ADGRA3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Adhesion G protein-coupled receptor A3 protein (101-381AA)
Target	ADGRA3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IWK6
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...