




Product Datasheet

CD38 antibody (orb34265)

Description	Rabbit polyclonal antibody to CD38
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CD38
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34328
Uniprot ID	P28907
NCBI	NP_001766.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HepG2 cells ...



Immunofluorescence analysis of human hep...