

## Product Datasheet

### NOTCH 1 peptide (orb343855)

<b>Description</b>	NOTCH 1 peptide
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	DOT
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Note</b>	For research use only
<b>Application notes</b>	Control peptide should be used at 1.0 µg per 1.0 µl of antiserum in per assay. Notch-1 (cleaved N-term) peptide has been tested in dot blot.
<b>Purity</b>	Greater than 95% specific peptide
<b>Source</b>	Human
<b>Uniprot ID</b>	<b>P46531</b>
<b>NCBI</b>	<b>CAG33502.1</b>
<b>Expiration Date</b>	12 months from date of receipt.