

Product Datasheet

Huntington (phospho-S421) Control peptide (orb343859)

Description

Huntington (phospho-S421) Control

Conjugation

Unconjugated

Tested

WB

Applications
Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Liquid (sterile filtered)

Concentration

1.0 mg/mL

Storage

Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Intended for use as a control in peptide competition assays to demonstrate the binding properties of anti-Huntington pS421 antibody to bind Huntington protein in the presence or absence of control peptide. Control peptide should be used at 1.0 µg per 1.0 µg of antibody per assay.

Purity

Greater than 95% specific peptide

Source

Human

Dilution Range

WB: 1.0 µg per 1.0 µg of antibody

Expiration Date

12 months from date of receipt.