

Product Datasheet

MAD2L1 peptide (orb343874)

Description	MAD2L1 peptide
Conjugation	Unconjugated
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Liquid (sterile filtered)
Concentration	1.0 mg/mL
Storage	Store vial at -20 °C prior to opening. Aliquot contents and freeze at -20 °C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4 °C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	MAD2L1 Control peptide should be used at 1.0 µg per 1.0 µl of antiserum in per assay.
Purity	Greater than 95% specific peptide
Source	Human
Expiration Date	12 months from date of receipt.