

Product Datasheet

H3 (phospho-Ser28) antibody (orb34390)

Description	Rabbit polyclonal antibody to H3 (phospho-Ser28). The protein belongs to the histone H3 family, and ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H3C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15404
Uniprot ID	P68431
NCBI	NP_066298.1, NP_003522.1, NP_003523.1, NP_003528.1, NP_003520.1, NP_003527.1, NP_003521.2, NP_003525.1, NP_003526.1, NP_003524.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

