

Product Datasheet

AKT (phospho-Ser129) antibody (orb34400)

Description	Rabbit polyclonal antibody to AKT (phospho-Ser129)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	AKT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55686
Uniprot ID	P31749
NCBI	NP_001014432.1 , NP_005154.2 , NP_001014431.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

