

Product Datasheet

C20orf196 antibody (orb34404)

Description	C20orf196 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein C20orf196 protein (1-205AA)
Target	C20orf196
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYI0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.