

Product Datasheet

APP (phospho-Ser730) antibody (orb34406)

Description	Rabbit polyclonal antibody to APP (phospho-Ser730). The protein belongs to the APP family, and is hi...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	APP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	86943
Uniprot ID	P05067
NCBI	NP_001191232.1, NP_958816.1, NP_001129601.1, NP_000475.1, NP_001129488.1, NP_001191231.1, NP_958817.1, NP_001129603.1, NP_001191230.1, NP_001129602.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

