



Product Datasheet

Caspase 6 (phospho-Ser257) antibody (orb34416)

Description	Rabbit polyclonal antibody to Caspase 6 (phospho-Ser257)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CASP6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33310
Uniprot ID	P55212
NCBI	NP_001217.2 , NP_116787.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in liver cell l...