

Product Datasheet

IZUMO1 antibody (FITC) (orb34427)

Description	IZUMO1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Izumo sperm-egg fusion protein 1 protein (22-292AA)
Target	IZUMO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYV9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.