

Product Datasheet

Protein A Agarose antibody (Agarose) (orb344346)

Description	Protein A Agarose antibody (Agarose)
Species/Host	Goat
Conjugation	Agarose
Immunogen	Protein A (Staphylococcus aureus.)
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Suitable for molecular biology, enzyme purification or to remove unwanted cross-reactivity.
Isotype	IgG
Clonality	Polyclonal
Purity	This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of Protein A.
Expiration Date	12 months from date of receipt.