

Product Datasheet

Avidin antibody (Biotin) (orb344992)

Description

Avidin antibody (Biotin)

Species/Host

Rabbit

Reactivity

Bacteria

Conjugation

Biotin

Tested

ELISA, IHC, WB

Applications

Immunogen

Avidin (Hen Egg White)

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Lyophilized

Concentration

1.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Avidin Antibody Biotin Conjugated has been tested by ELISA and western blot and is suitable for immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring lot-to-lot consistency.

Isotype

IgG

Clonality

Polyclonal

Purity

Anti-Avidin Antibody Biotin Conjugated is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Rabbit Serum, as well as purified and partially purified Avidin [Hen Egg White].

Uniprot ID


P02701

Dilution Range

ELISA: 1:2,000 - 1:10,000, IHC: 1:250 - 1,000, WB: 1:500 - 1:2,500

Expiration Date

12 months from date of receipt.



Western blot analysis of
Lane 1:
Avidin....