

Product Datasheet

SMAD4 (phospho-Thr277) antibody (orb34523)

Description Rabbit polyclonal to SMAD4 (phospho-Thr277)

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target SMAD4

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 60439

Uniprot ID [Q13485](#)


NCBI [NP_005350.1](#)

Dilution Range IF/ICC: 1:10-50, DOT: 1:500, IHC-P: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...