




Product Datasheet

STAT3 (phospho-Tyr705) antibody (orb34529)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	STAT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88068
Uniprot ID	P40763
NCBI	NP_998827.1, NP_644805.1, NP_003141.2
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formali...



Western blot analysis of STAT3 (phospho-...