


Product Datasheet

MAPKAPK5 antibody (orb34554)

Description	MAPKAPK5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human MAP kinase-activated protein kinase 5 protein (200-471AA)
Target	MAPKAPK5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IW41
Dilution Range	IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...