

Product Datasheet

JUN (phospho-Thr243) antibody (orb34555)

Description	Rabbit polyclonal to JUN (phospho-Thr243)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	JUN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35676
Uniprot ID	P05412
NCBI	NP_002219.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

