

Product Datasheet

CDH26 antibody (orb34576)

Description	CDH26 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Cadherin-like protein 26 protein (1-124AA)
Target	CDH26
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IXH8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.