

Product Datasheet

RFP antibody (orb345939)

Description

RFP antibody

Reactivity

Other

Conjugation

Unconjugated

Tested

WB

Applications
Preservatives

Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!

Form/Appearance

n/a

Storage

See kit insert for complete instructions.

Note

For research use only

Application notes

Use Biorbyt' Anti-RFP Chemiluminescent Kit for Western Blotting for detection of RFP-tagged recombinant proteins by western blot. This kit is useful for both western blotting and dot blotting methods. Please read the entire product insert prior to use.

Clonality

Polyclonal

Purity


Biorbyt' Chemiluminescent Western Blot Kit for RFP combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol chemiluminescent substrate for the detection of recombinant proteins containing wild type Red Fluorescent Protein (DsRed or RFP) and its variants such as mBanana, mCherry, DsRed2, E2-Crimson, mOrange, mOrange2, mPlum, mRaspberry, mStrawberry and tdTomato. The Chemiluminescent Western Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains sufficient substrate for up to 30 mini blots at 7.5 x 8 cm² (1,800 cm²) and is stable for at least 1 year when stored as indicated.

Uniprot ID
Q9U6Y8
Dilution Range

WB: User Optimized

Expiration Date

Please enquire.



Western blot analysis of HeLa cell-deriv...