

Product Datasheet

HRP Western Blot Mouse IgG antibody (Peroxidase) (orb345943)

Description	HRP Western Blot Mouse IgG antibody (Peroxidase)
Species/Host	Goat
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	1.0 mg/ml
Storage	See kit insert for complete instructions.
Note	For research use only
Application notes	Western Blot Analysis
Clonality	Polyclonal
Dilution Range	ELISA: 1:10,000-1:100,000, IHC: 1:500-1:5,000, WB: 1:5,000-1:40,000
Expiration Date	12 months from date of receipt.