



Product Datasheet

SLC38A6 antibody (orb34598)

Description	SLC38A6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Probable sodium-coupled neutral amino acid transporter 6 protein (132-169AA)
Target	SLC38A6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IZM9
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse kidney(la...



Immunofluorescence analysis of A431 cell...