

Product Datasheet

Bovine IgG Alkaline Phosphatase Antibody (orb346005)

Description

Bovine IgG Alkaline Phosphatase

Conjugation	AP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Liquid (sterile filtered)
Concentration	1.0 mg/mL
Storage	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.
Note	For research use only
Application notes	Bovine IgG whole molecule alkaline phosphatase conjugation can be utilized as a control or standard reagent in Western Blotting and ELISA experiments.
Isotype	IgG
Purity	This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by conjugation and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Bovine IgG, anti-Bovine Serum and anti-Alkaline Phosphatase (calf intestine).
Source	Bovine
Expiration Date	12 months from date of receipt.