



Product Datasheet

HSPB1 (phospho-Ser78) antibody (orb34603)

Description	Rabbit polyclonal to HSPB1 (phospho-Ser78)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	HSPB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22783
Uniprot ID	P04792
NCBI	NP_001531.1
Dilution Range	WB: 1:500, DOT: 1:500
Expiration Date	12 months from date of receipt.



WB analysis of Hela cell line using HSPB...