

Product Datasheet

Chicken IgG Agarose Antibody (orb346048)

Description	Chicken IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	IP, SDS-PAGE
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Chicken IgG Agarose has been tested by SDS-PAGE and is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product is an purified Chicken IgG isolated from normal Chicken Serum coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of anti-Chicken IgG.
Source	Chicken
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.