

Product Datasheet

Guinea Pig IgG Agarose Antibody (orb346132)

Description	Guinea Pig IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Guinea Pig IgG Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product was prepared by coupling purified Guinea Pig IgG antibody to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 10 mg of anti-Guinea Pig IgG antibody.
Source	Guinea Pig
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.