

Product Datasheet

Human IgG Agarose Antibody (orb346215)

Description	Human IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	SDS-PAGE
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	2.5 mg/mL
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Human IgG Agarose has been tested by SDS-Page and is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product was prepared by coupling purified Human IgG antibody to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of anti-Human IgG antibody.
Source	Human
Expiration Date	12 months from date of receipt.