

## Product Datasheet

### Rat IgG Agarose Antibody (orb346335)

<b>Description</b>	Rat IgG Agarose Antibody
<b>Conjugation</b>	Agarose
<b>Tested Applications</b>	IP
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Form/Appearance</b>	Suspension of agarose beads
<b>Storage</b>	Store vial at 4° C prior to opening. DO NOT FREEZE.
<b>Note</b>	For research use only
<b>Application notes</b>	Rat IgG Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
<b>Isotype</b>	IgG
<b>Purity</b>	This product is normal Rat IgG coupled to activated agarose. A single precipitin arc was observed against anti-Rat Serum when assayed by immunoelectrophoresis prior to coupling to the beads. Sufficient capacity is provided to bind a minimum of 10 mg of anti-Rat IgG.
<b>Source</b>	Rat
<b>Dilution Range</b>	IP: User Optimized
<b>Expiration Date</b>	12 months from date of receipt.