

Product Datasheet

PI4KAP1 antibody (HRP) (orb34637)

Description	PI4KAP1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Putative inactive phosphatidylinositol 4-kinase alpha-like protein P1 protein (1-262AA)
Target	PI4KAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N8J0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.