

Product Datasheet

Human AKT1 (phosphatase treated) protein (orb346472)

Description

Human AKT1 (phosphatase treated) protein

Conjugation

Unconjugated

Tested

ELISA, SDS-PAGE, WB

Applications

Preservatives

None

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Store vial at -70° C prior to use. Thaw only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. For long term storage we recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL) . For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature.

Note

For research use only

Application notes

Human AKT1 (phosphatase treated) recombinant protein has been tested in SDS-Page and western blot and is suitable as a control for polyclonal or monoclonal anti-AKT1 in immunological assays. Akt1 recombinant protein is dephosphorylated and in an inactive state. It is well suited as a control for anti-AKT pT308 which detects phosphorylated T308 residue, and anti-AKT pS473, that detects phosphorylated S473. For western blot use at 50 ng or less. For other assays concentration is user optimized.

Purity

Recombinant protein corresponds to amino acids 1 to 480 of mature human AKT1; Akt isoform 1. The recombinant protein contains a polyhistidine affinity tag at the amino terminus. Purity is greater than 90% as determined by reducing and non-reducing SDS-PAGE and by analytical HPLC.

Source

Human

Uniprot ID

P31749

NCBI


62241011

Dilution Range


ELISA: 50ng, WB: 50ng

Expiration Date

6 months from date of receipt.



Western blot analysis of Human Recombina...



Western blot analysis of Human Recombina...