



Product Datasheet

PAK1 (phospho-Thr423) antibody (orb34657)

Description	Rabbit polyclonal antibody to PAK1 (phospho-Thr423)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	PAK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60647
Uniprot ID	Q13153
NCBI	NP_002567.3 , NP_001122092.1
Dilution Range	IF/ICC: 1:10-50, DOT: 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...