

Product Datasheet

Chicken IgG antibody (Agarose) (orb346575)

Description	Chicken IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Gallus
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	10 mg/cc
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Chicken IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays.
Isotype	IgG
Clonality	Polyclonal
Purity	This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Chicken IgG.
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.