

Product Datasheet

Guinea Pig IgG antibody (Agarose) (orb346605)

Description	Guinea Pig IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Guinea pig
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Guinea Pig IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays.
Isotype	IgG
Clonality	Polyclonal
Purity	This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Guinea Pig IgG.
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.