

Product Datasheet

Golden Syrian Hamster IgG antibody (Agarose) (orb346615)

Description	Golden Syrian Hamster IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Golden Hamster
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Golden Syrian Hamster IgG Agarose Conjugate is useful for immunoprecipitation and other immunological assays.
Isotype	IgG
Clonality	Polyclonal
Purity	This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Hamster IgG.
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.