

Product Datasheet

Rat IgG antibody (Agarose) (orb346694)

Description	Rat IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Rat
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.1% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	4.4mg/mL
Storage	Store Anti-Rat IgG [H&L] (Goat) Antibody Agarose Conjugated vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Rat IgG [H&L] (Goat) Antibody Agarose Conjugated has been tested by SDDS-PAGE and is applicable to purification of Rat IgG from sample solution.
Isotype	IgG
Clonality	Polyclonal
Purity	Anti-Rat IgG [H&L] (Goat) Antibody Agarose Conjugated is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Rat IgG.
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.